# PROFI CARE®

## Bedienungsanleitung / Garantie

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## Instruction Manual

Thank you for choosing our product. We hope you will enjoy using the appliance.

### Symbols in these Instruction Manual

Important information for your safety is specially marked. It is essential to comply with these instructions in order to avoid accidents and prevent damage to the appliance:

#### WARNING:

This warns you of dangers to your health and indicates possible injury risks.

#### CAUTION:

This refers to possible hazards to the appliance or other objects.

**NOTE:** This highlights tips and information.

## Contents

Overview of the Components .....	3
General Notes .....	42
Special Safety Information Concerning the Rechargeable Built-in Battery .....	42
Toothbrush .....	43
Special Safety Instructions for this Appliance .....	43
Overview of the Components / Scope of Delivery .....	44
Charging the Appliance .....	44
Notes for Use .....	44
Usage .....	45
Cleaning .....	45
Oral Irrigator .....	45
Special Safety Instructions for this Appliance .....	45
Overview of the Components / Scope of Delivery .....	46
Charging the Appliance .....	46
Notes for Use .....	46

Initial Operation .....	47
Usage .....	47
Cleaning .....	47
Technical Data .....	47
Disposal .....	48

## General Notes

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this appliance to other people, please also pass on the operating instructions.

- The appliance is designed exclusively for private use and for the envisaged purpose. This appliance is not fit for commercial use.
- The charging station is suitable for use in dry rooms only.
- Do not touch the plugged-in charging station with wet hands!
- Turn off the oral irrigator and toothbrush when not in use, when attaching accessories, when cleaning or in case of malfunction.
- Do not use a damaged charging station.
- Do not repair the appliance yourself, but visit an authorized specialist.
- The appliance and the mains lead have to be checked regularly for signs of damage. If damage is found the appliance must not be used.
- Use only original spare parts.
- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

#### WARNING:

Do not allow small children to play with the foil.  
**There is a danger of suffocation!**

## Special Safety Information Concerning the Rechargeable Built-in Battery

### CAUTION:

- The oral irrigator and the toothbrush each contain a lithium-ion accumulator that must **not** be replaced by the user.
- Risk of explosion in case of improper replacement of the rechargeable battery. Replace only by the same or an equivalent type.

 **CAUTION:**

- Do not expose the rechargeable built-in battery to hard impacts or pressure. Do not drop the oral irrigator and the toothbrush!
- Do not expose the appliance to intense heat such as sunlight. Also keep the appliance away from heat sources like radiators, stoves, etc.!
- The ambient temperature should be between 0 °C and max. 35 °C. Heat damages the rechargeable battery. Increasingly lower temperatures reduce the capacity and the life of the rechargeable battery.
- Do not charge the rechargeable battery unsupervised.
- Disconnect the charging station from the power outlet if you are not going to use the oral irrigator and the toothbrush for a longer time.

## Toothbrush

### **Special Safety Instructions for this Appliance**

 **WARNING of electrical shock!**

- Keep the appliance dry!
  - Keep the charging station dry.
  - The appliance is intended for use with safety extra-low voltage. Connect only to a power supply according to the specifications on the rating label.
  - Do not replace any parts. Do not perform any modifications. Risk of electric shock!
- If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
- Do not use the appliance without a brush attachment.

- This appliance serves for the dental and oral hygiene of adult persons.
- This appliance can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.

### WARNING:

This appliance is not suited for children under the age of 3 as small parts could be broken off and swallowed.

- Children may not play with the appliance.
- Cleaning and user maintenance shall not be made by children without supervision.

### Overview of the Components / Scope of Delivery

#### Illustration A

- 1 Driving axle
- 2 POWER  $\phi$  button
- 3 Indicator light "C"
- 4 Indicator light "S"
- 5 Indicator light "M"
- 6 Indicator light "  "
- 7 Brush head
- 8 Brush attachment
- 9 Coloured marking ring
- 10 Slots for brush attachments
- 11 Charging surface for toothbrush
- 12 Charging station

#### No Illustration

- 2 replacement brushes

### Charging the Appliance

1. Connect the charging station to a correctly installed earthed socket.
2. Place the appliance on the charging station. The  control lamp turns on during charging
3. Charge the appliance for 24 hours before the first use.

#### NOTE:

- Put the toothbrush back to the charging station after use. Thus, it is always ready for use, even if you use it several times a day.
- Even when switched off, the rechargeable battery discharges continuously. If the appliance is out of use for a longer period of time, charge the rechargeable battery at least once every 4 to 8 weeks.

### Notes for Use

#### Brush Attachments

Delivery includes 3 brush attachments with different colour markings. These can be assigned to different users. 2 brush attachments can be stored in the slot of the charging station.

#### When the Gums are Bleeding ...

It is normal to experience slight gum bleeding when using the toothbrush for the first few times. Should your gums bleed continually, please consult your dentist.

#### $\phi$ Button

The  $\phi$  button turns the appliance on and off. You may also use the  $\phi$  button to select the programs.

Pressing the Button	Indicator Lamp	Function
Press 1x	C (Clean)	Turn appliance on
		For daily use
Press 2x	S (Sensitive)	Gently cleanses your teeth or other sensitive areas
Press 3x	M (Massage)	For a gentle gum stimulation
Press 4x	None	Turn appliance off

### Usage

- Place the brush attachment with the brush head facing the  $\Phi$  button. Press the brush attachment onto the appliance as far as it will go.
- Moisten the brush head and apply some toothpaste.
- Rinse your mouth with water.
- Lead the brush to the teeth before switching on the appliance via the  $\Phi$  button. You will thus prevent splatters.
- Place the brush against the edge of the gums. Brush the outsides, the insides and the occlusal surface of your teeth. Slowly guide the brush head from tooth to tooth. Only exert little pressure. In this way, clean each of the four quarters of your mouth for about 30 seconds. The appliance will indicate these intervals by a brief

interruption. After a total of 2 minutes, the recommended cleaning time by dentists, the toothbrush will turn off automatically.

#### **i** NOTE:

You may switch off the toothbrush at any time during the brushing process by using the  $\Phi$  button.

### Cleaning

#### **⚠** WARNING:

- Turn the appliance off. Remove the charging station plug from the socket when you clean the appliance.
- Do not submerge the appliance and the charging station in water. It could result in an electric shock or fire.

#### **⚠** CAUTION:

- Do not use any wire brush or other abrasive objects.
- Do not use any sharp or abrasive cleaners.
- You may wipe the housing of the appliance with a slightly damp cloth.
- Pull off the brush attachment. Clean it under running water after each use.
- Wipe the charging station with a slightly damp, not wet cloth.

## Oral Irrigator

### Special Safety Instructions for this Appliance

#### **⚠** WARNING of electrical shock!

- Keep the appliance dry!
  - Keep the charging station dry.
  - The appliance is intended for use with safety extra-low voltage. Connect only to a power supply according to the specifications on the rating label.
  - Do not replace any parts. Do not perform any modifications. Risk of electric shock!
- If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
  - Only use the appliance for the cleaning of the mouth.
  - Do not use the appliance without a jet nozzle.

- The appliance may only be used with cold or lukewarm water.
- Empty the water container after each application to prevent an accumulation of bacteria.
- When using braces, find out whether the usage of an oral irrigator is possible. We do not assume any liability for possible damages.
- In case of diseases in the mouth and throat area, you should ask your doctor whether the use of an oral irrigator is advisable.
- This appliance can be used by children over the age of 8 and by persons with reduced physical, sensory or mental capabilities or lack of experience and/or knowledge if they are supervised.
- Children may not play with the appliance.

### Overview of the Components / Scope of Delivery

#### Illustration B

- 1 Jet nozzle
- 2 Release button for the jet nozzle
- 3 ON/OFF button
- 4 MODE button
- 5 Indicator light "Normal"
- 6 Indicator light "Soft"
- 7 Indicator light "Pulse/CH"
- 8 Coloured marking ring
- 9 Suction hose
- 10 Filter
- 11 Sealing ring
- 12 Water tank
- 13 Sealing ring
- 14 Water tank cover
- 15 Charging surface for oral irrigator
- 16 Slots for jet nozzles
- 17 Charging station

#### No Illustration

- 2 replacement jet nozzles

### Charging the Appliance

1. Connect the charging station to a correctly installed earthed socket.
2. Place the appliance on the charging station. The indicator light "CH" (Charge) will blink during the charging procedure.
3. Charge the appliance for 24 hours before the first use.

#### **i** NOTE:

- If the appliance is fully charged, the indicator light "CH" will glow continuously.
- Charge the rechargeable battery if the indicator light "CH" blinks during or after using the appliance.
- Even when switched off, the rechargeable battery discharges continuously. If the appliance is out of use for a longer period of time, charge the rechargeable battery at least once every 4 to 8 weeks.

### Notes for Use

#### Jet Nozzles

3 jet nozzles with different colour markings are contained in the scope of delivery. These can be assigned to different users. 2 jet nozzles can be stored in the slot of the charging station.

#### When the Gums are Bleeding ...

It is normal for a normal bleeding to occur when using the oral irrigator the first few times. Should your gums bleed continually, please consult your dentist.

#### MODE Button

Select the desired water pressure setting using the MODE button:

**Normal:** Highest water pressure for a normal cleaning.

**Soft:** Lowest water pressure for gentle cleaning of sensitive gums.

**Pulse:** Water pressure at an interval for cleaning and massaging of the gums.

#### **i** NOTE:

- You can adjust the setting during use.

**i NOTE:**

- The higher the water pressure, the faster the water is used up.

### Initial Operation

1. Remove the oral irrigator from the charging station.
2. Attach a jet nozzle. Turn the jet nozzle so that the opening faces the front of the appliance.
3. Open the water tank cover on the back side of the appliance.
4. Fill the water tank with lukewarm water and close the water tank cover again. Make sure the sealing ring is fitted properly.
5. Select the "Normal" setting using the MODE button.
6. Turn on the appliance using the ON/OFF button.
7. Allow all water in the water tank to flow through the jet nozzle in order to remove any residue.
8. Turn off the appliance.
9. Briefly rinse off the jet nozzles under running water.

Your oral irrigator is now ready for use.

### Usage

Before using the water jet you should clean your teeth with a toothbrush.

1. Open the water tank cover on the back side of the appliance.
2. Fill the water tank with lukewarm water and close the water tank cover again.

**WARNING:**

Do not fill it with hot water!

3. In order to clean the jet nozzles, allow water to flow through the jet nozzle head (at maximum water pressure) into the wash basin before and after each time the appliance is used.

**NOTE:**

Do not hold the appliance at a too much slanting angle otherwise no water will be able to spray from the appliance.

4. Next, use the MODE button to select the desired setting "Normal", "Soft", "Pulse").
5. Lean over the wash basin. Guide the water jet along the gums, around the teeth. Hold the appliance upright as far as possible. Clean all four parts of your mouth in this way. Keep your mouth slightly open so the water can leak out.

**WARNING:**

Do not hold the jet nozzle too close to your gums or your teeth.

6. If the water becomes used up before you are finished with the cleaning, turn off the appliance using the ON/OFF button. Add more water.

7. To stop operation, turn the appliance off.

**NOTE:**

For security, the appliance automatically turns off after 2 minutes.

8. Please clean the appliance as described under "Cleaning".

### Cleaning

**WARNING:**

- Turn the appliance off. Remove the charging station plug from the socket when you clean the appliance.
- Do not submerge the appliance and the charging station in water. It could result in an electric shock or fire.

**CAUTION:**

- Do not use any wire brush or other abrasive objects.
- Do not use any sharp or abrasive cleaners.

1. Draw the water tank downwards.
2. If necessary, pour the remaining water out.
3. Clean the suction hose with a soft cloth.
4. Remove the filter and clean it carefully using warm running water.
5. Insert the filter in the hose once again.
6. Rinse out the water tank.
7. Guide the hose into the upper opening of the water tank.
8. Place the water tank in the appliance once again. Pay attention to the guides and make sure the sealing ring is fitted properly. The water tank must snap into place audibly.
9. Press the release button and remove the jet nozzle. Clean the jet nozzle under running warm water. Dry the jet nozzle.
10. Clean the appliance and the charging station with a damp cloth. Finally, wipe the parts until they are dry.

### Technical Data

Model:.....PC-DC 3031  
 Net weight: ..... approx. 622 g  
 Degree of protection against liquids: .....IPX7\*

#### Charging station

Protection class:..... II  
 Power supply:..... 100-240 V~, 50-60 Hz  
 Power consumption: ..... 1.2-2.8 W

**Toothbrush accumulator:** ..... 1 × 3.7 V, 650mAh Li-Ion  
 Running time: ..... 4.5 hours max.  
 Charging time: ..... approx. 24 hours

**Oral irrigator accumulator:** ..... 1 × 3.7 V, 1500mAh Li-Ion  
 Running time: ..... 60 minutes max.  
 Charging time: ..... approx. 24 hours

\* This appliance complies with IPX7. This means that the appliance is protected against temporary submersion.

The right to make technical and design modifications in the course of continuous product development remains reserved.

This appliance has been tested according to all relevant current CE guidelines, such as electromagnetic compatibility and low voltage directives, and has been constructed in accordance with the latest safety regulations.

## Disposal

### Dispose of Batteries/Rechargeable Batteries

The oral irrigator and the toothbrush each contain a lithium-ion accumulator. Used batteries/accumulators (rechargeable batteries) do not belong in household waste.

The battery must be removed before the product is recycled. Contact qualified personnel.

Information, where batteries/accumulators are to be disposed of can be obtained from your local or municipal authorities.

### Meaning of the "Dustbin" Symbol

Protect our environment: do not dispose of electrical equipment in the domestic waste.

Please return any electrical equipment that you will no longer use to the collection points provided for their disposal.

This helps avoid the potential effects of incorrect disposal on the environment and human health.

This will contribute to the recycling and other forms of re-utilisation of electrical and electronic equipment.

Information concerning where the equipment can be disposed of can be obtained from your local authority.

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